

# Merckodia

# Adiponectin ELISA

Directions for Use

**10-1193-01**

**REAGENTS FOR 96 DETERMINATIONS**

For Research Use Only  
Not for Use in Diagnostic Procedures





Manufactured by

**Merckodia AB**

Sylveniusgatan 8A  
SE-754 50 Uppsala  
Sweden

Merckodia 

## EXPLANATION OF SYMBOLS USED ON LABELS

 $\Sigma = 96$	Reagents for 96 determinations
	Expiry date
	Store between 2-8°C
	Lot No.

## INTENDED USE

Mercodia Adiponectin ELISA provides a method for the quantitative determination of human adiponectin in serum or plasma.

## SUMMARY AND EXPLANATION OF THE TEST

Adiponectin is also called: Acrp30 (30 kDa adipocyte complement-related protein), GBP28 (Gelatin-binding protein), adipoQ, apM1 (Adipose most abundant gene transcript 1) [1, 2].

Adiponectin is an adipocyte-secreted hormone, consisting of 244 amino acids with a molecular weight of approximately 30kDa (28-30kDa). It is one of the most abundant proteins in human blood, with circulating concentrations of 0.5-30 µg/ml, which accounts for approximately 0,01% of total plasma protein [2].

The protein consists of four domains: one globular C-terminal, one collagen-like N-terminal, one signaling peptide and one hyper variable domain. The globular domain has significant sequence and structural similarities to the complement factor C1q [2,3]. The globular domain also has structural similarities to TNF-α [3-5].

Adiponectin concentration is reversely associated with type 2 diabetes, coronary artery disease and obesity, all together called the metabolic syndrome. Adiponectin decreases blood glucose and free fatty acid serum concentrations and increases insulin sensitivity [7]. Adiponectin has also been shown to have anti-inflammatory effects [2].

Adiponectin has been suggested to exist in different forms in circulation: monomers, isolated globular forms (the globular domains), trimers, hexamers and larger oligomers [8-11]. Monomers are believed to associate in circulation to trimers through the globular domains. Trimers are associated to larger oligomers through the collagen-like domain [7].

However, recent studies indicate that adiponectin may not be present in circulation as monomers or isolated globular forms, but rather in multimeric structures. The studies have shown that the dominant forms of adiponectin that circulates in human blood are hexamers (LMW) and larger oligomers (HMW) [6, 12-14]. The LMW adiponectin levels does not seem to differ between insulin sensitive- and insulin resistant subjects, nor does LMW adiponectin differ between men and women. The increased levels of total adiponectin in insulin sensitive subjects and women were caused by increased amounts of HMW adiponectin. Both total and HMW adiponectin showed significant differences between the insulin sensitive- and insulin resistant subjects according to Lara-Castro et al. 2006 [6].

Several isoforms of adiponectin do circulate in blood. It is yet to be determined whether all isoforms are secreted by the adipocytes, whether there is a posttranscriptional assembly of HMW adiponectin in blood or whether the HMW form is secreted and degraded in blood. The individual metabolic significance of each adiponectin isoform also remains unclear [6].

## **PRINCIPLE OF THE PROCEDURE**

Mercodia Adiponectin ELISA is a solid phase two-site enzyme immunoassay. It is based on the sandwich technique in which two monoclonal antibodies are directed against separate antigenic determinants on the adiponectin molecule. During incubation, adiponectin in the sample react with anti-adiponectin antibodies bound to microtitration well. After washing, peroxidase conjugated anti-adiponectin antibodies are added and after the second incubation and a simple washing step that removes unbounded enzyme labelled antibody, the bound conjugate is detected by reaction with 3,3',5,5'-tetramethylbenzidine (TMB). The reaction is stopped by adding acid to give a colorimetric endpoint that is read spectrophotometrically.

## **WARNINGS AND PRECAUTIONS**

- For research use only. Not for use in diagnostic procedures. Not for internal or external use in humans or animals.
- The content of this kit and their residues must not be allowed to come into contact with ruminating animal or swine.
- The Stop Solution in this kit contains 0.5 M H<sub>2</sub>SO<sub>4</sub>. Follow routine precautions for handling hazardous chemicals.
- All patient samples should be handled as capable of transmitting infections.

## **MATERIAL REQUIRED BUT NOT PROVIDED**

- Pipettes for 20, 25, 50, 100, 200 and 1000 µl (repeat pipettes preferred for addition of enzyme conjugate solution, Substrate TMB and Stop Solution)
- Beakers and cylinders for reagent preparation
- Redistilled water
- Microplate reader (450 nm filter)
- Plate shaker (The recommended velocity is 700-900 cycles per minute, orbital movement)
- Microplate washing device

## REAGENTS

Each Mercodia Adiponectin ELISA kit (10-1193-01) contains reagents for 96 wells, sufficient for 42 samples and one calibrator curve in duplicate. For larger series of assays, use pooled reagents from packages bearing identical lot numbers. The expiry date for the complete kit is stated on the outer label. The recommended storage temperature is 2-8°C.

<b>Coated Plate</b> (Mouse monoclonal anti-human adiponectin)	1 plate 8-well strips	96 wells	Ready for use
For unused microtitration strips, reseal the bag using adhesive tape and store at 2-8°C for two months.			
<b>Calibrators 1, 2, 3, 4, 5</b> (Recombinant human adiponectin) Concentration stated on vial label.	5 vials	1000 µl	Ready for Use
<b>Calibrator 0</b> Color coded yellow	1 vial	5 ml	Ready for use
<b>Assay Buffer</b> Color coded red	1 vial	12 ml	Ready for use
<b>Sample Buffer 2X</b> Dilute with 50 ml redistilled water to make sample buffer. Color coded yellow Storage after dilution: 2-8°C for two months	1 bottle	50 ml	
<b>Enzyme Conjugate 11X</b> (Peroxidase conjugated mouse monoclonal anti-human adiponectin) <i>Note! Light sensitive!</i>	1 vial	1.3 ml	Preparation, see below
<b>Enzyme Conjugate Buffer</b> Color coded blue	1 vial	13 ml	Ready for use
<b>Wash Buffer 21X</b> Dilute with 800 ml redistilled water to make wash buffer. Storage after dilution: 2-8°C for two months	1 bottle	40 ml	
<b>Substrate TMB</b> (TMB) Colorless solution <i>Note! Light sensitive!</i>	1 bottle	22 ml	Ready for use
<b>Stop Solution</b> 0.5 M H <sub>2</sub> SO <sub>4</sub>	1 vial	7 ml	Ready for use

## Preparation of enzyme conjugate solution

Prepare the needed volume of enzyme conjugate solution by dilution of Enzyme Conjugate 11X, (1+10) in Enzyme Conjugate Buffer according to the table below. Mix gently.

When preparing enzyme conjugate solution for the whole plate or if the reagents are to be used within two months, pour all of the Enzyme Conjugate Buffer into the Enzyme Conjugate 11X vial.

Number of strips	Enzyme Conjugate 11X	Enzyme Conjugate Buffer
12 strips	1 vial	1 vial
8 strips	700 µl	7 ml
6 strips	500 µl	5 ml
4 strips	400 µl	4 ml

Storage after dilution: 2-8°C for two months.

## SPECIMEN COLLECTION AND HANDLING

### Serum

Collect blood by venipuncture, allow to clot. Separate the serum by centrifugation at 4 300 g for 15 minutes at 2-8°C. Specimen can be stored at 2-8°C up to 14 days. For longer periods, store samples at -20°C. Avoid repeated freezing and thawing.

### Plasma

Collect blood by venipuncture into tubes containing heparin, citrate or EDTA as anticoagulant, and separate the plasma fraction. Samples can be stored at 2-8°C up to 14 days. For longer periods store samples at -20°C. Avoid repeated freezing and thawing.

## PREPARATION OF SAMPLES

Samples should be diluted 1/101 v/v with sample buffer (20 µl sample + 2,0 ml sample buffer). Diluted samples can be stored at 2-8°C up to 14 days. *Note!* Buffers containing sodium azide (NaN<sub>3</sub>) can not be used for sample dilution.

## TEST PROCEDURE

All reagents and samples must be brought to room temperature before use. Prepare a calibrator curve for each assay run.

1. Prepare enzyme conjugate solution (according to the table on previous page), sample buffer, wash buffer and samples.
2. Prepare sufficient microplate wells to accommodate Calibrators and samples in duplicate.
3. Pipette 25  $\mu$ l each of Calibrators and samples into appropriate wells.
4. Add 100  $\mu$ l of Assay Buffer into each well.
5. Incubate on a plate shaker (700-900 rpm) for 1 hour at room temperature (18-25°C).
6. Wash plate 6 times with automatic plate washer  
or  
aspirate the reaction volume completely and fill each well with 350  $\mu$ l wash buffer. Aspirate liquid completely. Repeat 5 times. After final wash, invert and tap the plate firmly against absorbent paper.
7. Add 100  $\mu$ l of enzyme conjugate solution into each well.
8. Incubate on a plate shaker for 1 hour at room temperature (18-25°C).
9. Wash plate 6 times with automatic plate washer  
or  
aspirate the reaction volume completely and fill each well with 350  $\mu$ l wash buffer. Aspirate liquid completely. Repeat 5 times. After final wash, invert and tap the plate firmly against absorbent paper.
10. Add 200  $\mu$ l Substrate TMB into each well.
11. Incubate for 15 minutes at room temperature (18-25°C).
12. Add 50  $\mu$ l Stop Solution to each well.  
Place the plate on the shaker for approximately 5 seconds to ensure mixing.
13. Read optical density at 450 nm and calculate results.  
Read within 30 minutes.

*Note!* To prevent contamination between the conjugate and substrate, separate pipettes are recommended.

## INTERNAL QUALITY CONTROL

Internal serum pools with low, intermediate and high adiponectin concentrations should routinely be assayed as unknowns, and results charted from day to day. It is good laboratory practice to record the following data for each assay: kit lot number; reconstitution dates of kit components; OD values for the blank (sample buffer), Calibrators and controls.

## CALCULATION OF RESULTS

### Computerized calculation

The concentration of adiponectin is obtained by computerized data reduction of the absorbance for the Calibrators 1-5 versus the concentration using cubic spline regression.

### Manual calculation

1. Plot the absorbance values obtained for the Calibrators 1-5 against the adiponectin concentration on a log log paper and construct a calibrator curve.
2. Read the concentration of the unknown samples from the calibrator curve.
3. Multiply the concentration with the dilution factor.


### Example of results

Wells	Identity	A <sub>450</sub>	Mean conc. ng/ml	x101 µg/ml
1A-B	Calibrator 0	0.059/0.056		
1C-D	Calibrator 1 (5 ng/ml)*	0.106/0.102		
1E-F	Calibrator 2 (15 ng/ml)*	0.207/0.216		
1G-H	Calibrator 3 (50 ng/ml)*	0.563/0.559		
2A-B	Calibrator 4 (150 ng/ml)*	1.477/1.567		
2C-D	Calibrator 5 (300 ng/ml)*	2.602/2.681		
2E-F	Unknown 1	0.374/0.367	29.825	3.012
2G-H	Unknown 2	0.754/0.758	67.178	6.785
3A-B	Unknown 3	1.385/1.373	133.340	13.467

\*Exact concentration indicated on vial label.

### Example of calibrator curve

A typical calibrator curve is shown here. Do not use this curve to determine actual assay results.



## LIMITATIONS OF THE PROCEDURE

As with all diagnostic tests, a definitive diagnosis should not be based on the results of a single test, but should be made by a physician after all clinical findings have been evaluated. Grossly lipemic, icteric or haemolyzed samples do not interfere in the assay.

## EXPECTED VALUES

Good practice dictates that each laboratory establishes its own expected range of values.


Mercodia Adiponectin ELISA detects LMW(Hexamer 230kDa) and HMW(Oligomer >420 kDa) adiponectin, as determined by size exclusion gel chromatography.

The different multimeric forms of endogenous adiponectin were studied and separated in serum from a healthy individual by a three-step method; ammonium sulphate precipitation followed by ion exchange and gel filtration chromatography.

With ion exchange chromatography, proteins binds to the matrix with electrostatic forces causing separation since different proteins/isoforms have different total net charges or isoelectric points. The ion exchange column used was Mono Q 10/100GL (GE Healthcare). Triethanolamine buffer was used for eluting the proteins.

The isoforms of adiponectin has been shown to have different isoelectric points and post-translational patterns [15]. Proline hydroxylation and lysine hydroxylation/glycosylation are believed to have great importance on the assembly of the oligomers [12,15].


Five clearly distinct peaks were visible when the serum adiponectin was separated by ion exchange chromatography, indicating that the serum adiponectin analyzed presents at least five different post-translational patterns, yielding different isoelectric points, see figure 1 below. Pools A-E were further analyzed by size exclusion gel chromatography to determine the apparent size of the multimeric adiponectin forms.



**Figure 1** Elution profile of serum adiponectin from ion exchange chromatography and as identified by the Mercodia Adiponectin ELISA. Peak fractions were pooled to A, B, C, D and E.

Size exclusion gel filtration chromatography separates proteins according to apparent globular size. The gel filtration column that was used was HiLoad 16/60 Superdex 200 prep grade (GE Health care). PBS was used for eluting the proteins.

Three dominant multimetric forms were visible when the serum adiponectin was separated by size exclusion gel chromatography, with apparent sizes of 230 kDa, 420 kDa and > 600 kDa respectively, and interpreted as LMW (hexamer 230 kDa) and HMW (420 kDa and >600 kDa), see figure 2 below.



**Figure 2** Elution profile of serum adiponectin from size exclusion gel filtration chromatography and as identified by the Mercodia Adiponectin ELISA. Each pool (A,B, C, D and E) was analyzed separately and together as one pool (Pool).

In conclusion, the serum adiponectin analyzed displayed three dominant multimetric forms based on size, and five different forms based on isoelectric points, or total net charges.

## PERFORMANCE CHARACTERISTICS

### Detection limit

Detection limit is defined as the Capability of Detection according to ISO11843-Part 1. Capability of Detection should be seen as part of a method validation, rather than the lowest concentration that can be measured.

The detection limit is 1.25 (ng/ml) as determined by the methodology described in ISO11843-Part 4.

Concentration of samples with absorbance below Calibrator 1 should not be calculated, instead expressed as less or equal to ( $\leq$ ) the concentration indicated on the vial for Calibrator 1.

### Recovery

Recovery upon addition is 92-109% (mean 101%).

Recovery upon dilution is 89-111% (mean 98%).

### Hook effect

There is no existing hook effect.

### Precision

Each sample was analyzed in 4-replicates on 39 different occasions.

Sample	Mean value (ng/ml)	Coefficient of variation		
		within assay %	between assay %	total assay %
1	29.7	3.0	5.3	5.5
2	65.9	2.7	5.0	5.2
3	13.0	3.0	5.8	6.0

### Specificity

The following crossreactions have been found:

C1q	$\leq 0.007\%$
TNF- $\alpha$	n.s

## CALIBRATION

Adiponectin ELISA is calibrated against a highly purified, fully validated, commercial adiponectin preparation. The concentration of adiponectin is expressed in ng/ml.

## WARRANTY

The performance data presented here was obtained using the procedure indicated. Any change or modification in the procedure not recommended by Mercodia AB may affect the results, in which event Mercodia AB disclaims all warranties expressed, implied or statutory, including the implied warranty of merchantability and fitness for use. Mercodia AB and its authorized distributors, in such event, shall not be liable for damages indirect of consequential.

## REFERENCES

- [1] PDB, Protein Data Bank: [www.expasy.org](http://www.expasy.org), adiponectin: Q 15848 (21 oct.2005).
- [2] Meier U and Gressner AM (2004) Endocrine Regulation of Energy Metabolism: Review of Pathobiochemical and Clinical Aspects of Leptin, Ghrelin, Adiponectin, and Resistin. *Clin Chem* 50:1511-1525
- [3] Tsao TS, Lodish HF and Fruebis J (2002) ACRP30, a new hormone controlling fat and glucose metabolism. *Eur J Pharmacol* 440:213-221
- [4] Shapiro L and Scherer PE (1998) The crystal structure of a complement-1q family protein suggests an evolutionary link to tumor necrosis factor. *Curr Biol* 12:335-338
- [5] Bobbert T, Rochlitz H, Wegewitz U, Akpulat S, Mai K, Weickert MO, Mohlig M, Pfeiffer AF and Spranger J (2005) Changes of adiponectin oligomer composition by moderate weight reduction. *Diabetes* 54:2712-2719
- [6] Lara-Castro C, Luo N, Wallace P, Klein RL and Garvey T (2006) Adiponectin Multimeric Complexes and the Metabolic Syndrome Trait Cluster. *Diabetes* 55: 249-259
- [7] Duntas LH, Popovic V and Panotopoulos G (2004) Adiponectin: Novelities in Metabolism and Hormonal Regulation. *Nutr Neurosci* 7:195-200
- [8] Tsao TS, Lodish HF and Fruebis J (2002) ACRP30, a new hormone controlling fat and glucose metabolism. *J Biol Chem* 50:50810-50817
- [9] Nakashima R, Kamei N, Yamane K, Nakanishi S, Nakashima A and Kohno N (2006) Decreased Total and High Molecular Weight Adiponectin Are Independent Risk Factors for the Development of Type 2 Diabetes in Japanese-Americans. *Clin Endocrinol Metab* 91:3873-3877
- [10] Hara K, Horikoshi M, Yamauchi T, Yago H, Miyazaki O, Ebinuma H, Imai Y, Nagai R, and Kadowaki T (2006) Measurement of the High-Molecular Weight Form of Adiponectin in Plasma Is Useful for the Prediction of Insulin Resistance and Metabolic Syndrome. *Diabetes Care* 29:1357-1362

- [11] Nakano Y, Tobe T, Choi-Miura NH, Mazda T and Tomita M (1996) Isolation and Characterization of GBP28, a Novel Gelatin-Binding Protein Purified from Human Plasma. *J Biol Chem* 120:803-812
- [12] Bodles AM, Banga A, Rasouli N, Ono F, Kern PA and Owens RJ (2006) Pioglitazone increases secretion of high-molecular-weight adiponectin from adipocytes. *Am J Physiol Endocrinol Metab* 291:E1100-E1105
- [13] Halperin F, Beckman JA, Patti ME, Trujillo ME, Garvin M, Creager MA, Scherer PE and Goldfine AB (2005) The role of total and high-molecular-weight complex of adiponectin in vascular function in offspring whose parents both had type 2 diabetes. *Diabetologia* 48:2147-2154
- [14] Fisher FF, Trujillo ME, Hanif W, Barnett AH, McTernan PG, Scherer PE, and Kumar S (2005) Serum high molecular weight complex of adiponectin correlates better with glucose tolerance than total serum adiponectin in Indo-Asian males. *Diabetologia* 48:1084-1087
- [15] Wang Y, Lu G, Wong WP, Vliegenthart JF, Gerwig GJ, Lam KS, Cooper GJ and Xu A (2004) Proteomic and functional characterization of endogenous adiponectin purified from fetal bovine serum. *Proteomics* 4:3933-3942





## SUMMARY OF PROTOCOL SHEET

Add Calibrators and samples	25 $\mu$ l
Add Assay Buffer	100 $\mu$ l
Incubate	1 hour at 18-25°C on a plate shaker
Wash plate with wash buffer	6 times
Add enzyme conjugate solution	100 $\mu$ l
Incubate	1 hour at 18-25°C on a plate shaker
Wash plate with wash buffer	6 times
Add Substrate TMB	200 $\mu$ l
Incubate	15 minutes at 18-25°C
Add Stop Solution	50 $\mu$ l Shake for 5 seconds to ensure mixing
Measure A450	Evaluate results